

# Hydrall AA

# Aplicación 047




innovación  
tecnológica  
para  
laboratorio

# Rafer

## Determinación de mercurio en orina mediante CVAAS

### Introduction

Mercury levels in urine can be used to help diagnose recent mercury exposure and to evaluate patient response to chelation therapy. This method was developed to minimize sample pretreatment without loss in accuracy and utilizes bromine monochloride as an oxidant.



### Reagents

**Bromine monochloride (BrCl)**—In a fume hood, dissolve 27 g of reagent grade KBr in 2.5 L of low-Hg HCl. Place a clean magnetic stir bar in the bottle and stir for

approximately 1 h in the fume hood. Slowly add 38 g reagent grade  $\text{KBrO}_3$  to the acid while stirring. When all of the  $\text{KBrO}_3$  has been added, the solution color should change from yellow to red to orange. Loosely cap the bottle, and allow to stir another hour before tightening the lid.

**Hydroxylamine hydrochloride**—Dissolve 300 g of  $\text{NH}_2\text{OH}\cdot\text{HCl}$  in reagent water and bring to 1.0 L. This solution may be purified by the addition of 1.0 ml of  $\text{SnCl}_2$  solution and purging overnight at 500 ml/min with Hg-free  $\text{N}_2$ .

**Stannous chloride** – Add 100 ml of concentrated (12N) hydrochloric acid to about 500 ml of de-ionized water. Next, add 20 grams of anhydrous stannous chloride and swirl to dissolve. Add de-ionized water to 1.0 L.

**Rinse solution** – Add 20 ml of concentrated (12N) hydrochloric acid to about 500 ml of deionized water. Next, add de-ionized water to 1.0L.

## Sample Pretreatment

Dispense 0.5 ml of urine sample into 15 ml polypropylene test tubes. Add to each sample 9.9 ml of de-ionized water (18 MOhm) and 0.1 ml of bromine monochloride (BrCl) solution. Let mixture stand for two hours and immediately before analysis add 0.010 ml 30% hydroxylamine hydrochloride.

## Analysis


Calibrate the spectrometer with aqueous standards in the range of 0-10 ug/L. We used standards of 0, 0.2, 0.5, 1.0, 5.0, and 10.0 ug/L. Table I shows the operational parameters for the Hydra AA.

Parameter	Value
Carrier Gas	0.3 LPM
Pump Speed	7 ml/min
Rinse Time	60 sec.
Uptake Time	20 sec.
Integration Time	45 sec.
Fit Type	Wt. Linear
Method	CVAFS

*Table I: Operational Parameters*

## Results

Figure I shows the calibration that was obtained using aqueous mercury standards stabilized with 2% hydrochloric acid (same as rinse solution). The correlation coefficient (Rho) for the curve using weighted linear fit was 0.99964 and replicate readings exhibited a precision typically better than 1%.



*Figure 1: Calibration*

Reference urine samples were obtained from Bio-Rad Laboratories, Irvine, CA. Mercury levels for each sample appear below in Table II together with the values obtained with this method.

Sample	Lot	Mean	Range	Measured
Lypchocek 1	69091	47 ug/L	38-57ug/L	49.8 ug/L
Lypchocek 2	69092	147	117-176	154

**Table II: Reference Samples**

## Conclusions

The Hydra AA has the sensitivity and precision to determine mercury at the typical levels found in healthy individuals (10-20ug/L) as well as the accuracy to determine mercury at elevated levels found in patients recently exposed to mercury or on chelation therapy.

Bromine monochloride digestion effectively releases the mercury in urine for determination with calculated recoveries of 106% (Lypchek1) and 105% (Lypchek2) without the need for heating or extended reaction times.