

## Determination of bromine index of aromatic hydrocarbon

Sheet **GT200-OF033E** PetroleumMethod **Polarization titration**Automatic Titrator Model GT-200  
(GT0EF)Related standard ASTM 5776-07  
Standard Test Method for  
Bromine Index of Aromatic  
Hydrocarbons by  
Electrometric Titration

Electrodes:

Apparatus Double platinum electrode (GTRE10B)

\*Optional unit

1) Analog pack PS board for polarization  
and conductivity (GTEPSK)

2) Cooling jacket for tall beaker

Titration mode INF, Detection: mV, Differential  
polarographic titration/IPOL, Applied  
voltage: 300mV

※ \*This application sheet is provided as reference, and does not assure the measurement results. Please consider analysis environment, external factors and sample nature for optimal conditions before the measurement.

### Outline

A bromine number is a number which expresses the content of unsaturated components contained in petroleum or the like by the number of grams of bromine added to unsaturated components in a 100g sample (gBr<sub>2</sub>/100g). A bromine index is an index which expresses 1 gBr<sub>2</sub>/100g of a bromine number as a bromine index 1000, which corresponds to the number of milligrams of bromine added to unsaturated components.

### Reagents

[Titrant] ■ 0.05mol/L (N/10)-bromine solution (Volumetric analysis grade/Commercial item)  
[Titration solvent] Mixture of 714ml of acetic acid (Special grade), 134ml of 1-Methyl-2-pyrrolidinone, 134ml of methanol and 18ml of sulfuric acid (Sulfuric acid:Water = 1:5)

### Analytical procedure

[Blank measurement]

- (1) Collect 150ml titration solvent using a measuring cylinder and add it into a 200ml tall beaker.
- (2) Set the above-mentioned beaker in a cooling jacket and titrate with 0.05mol/L (N/10)-bromine solution while cooling the beaker down to 0-5 degrees C using a cold water circulation device or the like.

[Sample titration]

- (1) Collect 30ml titration solvent using a measuring cylinder and add it into a 200ml tall beaker.
- (2) Set the above-mentioned beaker on a balance and weigh out 20g of the sample.
- (3) Collect 120ml titration solvent using a measuring cylinder and add it into the beaker.
- (4) Set the above-mentioned beaker in a cooling jacket and titrate with 0.05mol/L (N/10)-bromine solution while cooling the beaker down to 0-5 degrees C using a cold water circulation device or the like.

[Calculation]

$$\text{Bromine index} = (A1 - BL) \times K1 \times 7990/S$$

- A1: Titration volume of 0.05mol/L-bromine solution at sample titration (ml)  
 BL: Titration volume of 0.05mol/L-bromine solution at blank measurement (ml)  
 K1: Normality of bromine solution (0.1)  
 7990: Atomic weight of bromine (79.90) × 100. (Converted into the value per 100g)  
 S: Sample volume (g)

### Other requirement

- Make sure to confirm labels and safety data sheets of reagents and gases used for the measurement and handle them with enough care.
- Wear protective equipment (eye protector, gloves and others) when handling reagents.
- Refer to the instruction manual of the analog pack polystyrene substrate for the setting method of applied current.
- When using a cooling jacket, put approximately 40ml of water in it as a medium. As titration is performed while cooling the water down to 0-5 degrees C, prepare a cold water circulation device or the like separately.

### Measurement results

	Sample size (g)	Titration volume (ml)	Bromine index
1	20.2574	5.3614	210.6
2	20.0714	5.2426	207.8
3	20.5043	5.3615	208.1


Blank: 0.0218ml

Nos. of data	(n)	3
Average		208.8
Standard deviation	(SD)	1.54
Relative standard deviation	(RSD%)	0.74

Bromine indexes of aromatic hydrocarbon were measured using GT-200. The average of three measurements was 208.8 and the relative standard deviation (RSD %) was 0.74%. GT-200 can measure bromine indexes with good repeatability.

Measurement : 2015/05/08 14:06  
 Sample name : Titration solvent 150ml

Type : Sample Titr  
 Sample size (S) : 150 [ml]



C1 : 0.0218 [ml]

A1 : 0.0218 [ml] 82 [mV]

Initial potential (Pi) : 19 [mV]  
 Start : 0 [ml] 19 [mV]  
 End : 0.21 [ml] 355 [mV] Time : 2'39"


Run file No. : 5 Bromine index

Titration filr No. : 9 Blank of Bromine index \*Run file and Titration file parameters are set for each analysis item

Mode : INF End1, End1 Width : 100 [mV] ± 100 [mV]  
 Detect : mV(P)  
 BRT No. : 1  
 Reagent : 35  
 WTint : 10 [sec]  
 Vup : 10 [μl]  
 Vlow : 10 [μl]  
 dE : 3 [mV]  
 dT : 5 [sec]  
 DL : 500 [mV/ml]  
 DetCnt : 3  
 Vmax : 5 [ml]  
 Vover : 0.1 [ml] C1 : A1 [ml]

Reagent name (Reag) : Bromine solution Equivalent (E) : 2 Molarity (M) : 0.05 [Mol/l]  
 Factor (f) : 1.005

Buret Injection Speed : 500 [ul/sec]



Run file No. : 5 Bromine index  
 Titration file No. : 5 Bromine index      \*Run file and Titration file parameters are set for each analysis item  
 Mode : INF      End1, End1 Width : 500 [mV] ± 500 [mV]  
 Detect : mV(P)  
 BRT No. : 1  
 Reagent : 35  
 WTint : 10 [sec]  
 Vup : 200 [μl]  
 Vlow : 200 [μl]  
 dE : 3 [mV]  
 dT : 5 [sec]  
 DL : 500 [mV/ml]  
 DetCnt : 6      C1 : (A1-BL)\*K1\*7990/S  
 Vmax : 100 [ml]  
 Vover : 1 [ml]

Reagent name (Reag) : Bromine solution      Equivalent (E) : 2      Molarity (M) : 0.05 [Mol/l]  
 Factor (f) : 1.005      Blank (BL) : 0.0218 [ml]      Coefficient1(K1) : 0.1

Buret Injection Speed : 500 [ul/sec]