INSTRUCTIONS FOR USE

INTENDED USE

CD-Chex Plus® BC is a stabilized preparation of normal human peripheral blood intended for use as a quality control material for Beckman Coulter® flow cytometry instruments.

SUMMARY/PRINCIPLES

Immunophenotyping by flow cytometry provides a rapid and accurate assessment of the frequency and type of leukocytes present in a blood sample. Quantitation of peripheral blood leukocyte subsets has been shown to be useful in the evaluation of immunodeficiency, autoimmune, and oncologic disorders. CD-Chex Plus BC control cells possess surface antigens detectable with monoclonal antibodies. When these cells are stained with fluorescent antibodies and analyzed by flow cytometry, they provide a reference for normal peripheral blood leukocytes.

Subsets of leukocytes can be distinguished on the basis of light scatter properties as well as the presence of cell surface antigens. CD-Chex Plus BC is designed to represent normal peripheral blood leukocytes where lymphocytes, monocytes, and granulocytes are distinguishable on the basis of light scatter properties.

REAGENTS

This product contains stabilized human red blood cells and white blood cells in a preservative medium.

PRECAUTIONS

- For In Vitro Diagnostic Use Only.
- CAUTION:** All blood products should be treated as potentially infectious. Source material from which this product was derived was found negative when tested in accordance with current FDA required tests. No known test methods can offer assurance that products derived from human blood will not transmit infectious agents. See the Instructions (IFU) tab under Resources on the product page at www.streck.com for specific FDA required blood tests.
- This product should not be disposed in general waste, but should be disposed with infectious medical waste. Disposal by incineration is recommended.
- This product is intended for use as supplied. Adulteration by dilution or addition of any materials to the product vial invalidates any diagnostic use of the product.
- CD-Chex Plus BC should not be used as a calibrator.

STORAGE AND STABILITY

CD-Chex Plus BC cells are stable through the expiration date when stored at 2 °C to 10 °C. After opening, CD-Chex Plus BC is stable throughout the open-vial dating, as indicated on the assay sheet, when stored at 2 °C to 10 °C. **DO NOT FREEZE.**

INDICATIONS OF PRODUCT DETERIORATION

If CD-Chex Plus BC values are not within the expected range on the assay sheet or granulocytes show a loss of forward scatter (FSC):

- Review the operating procedure of the instrument.
- Assay an unopened vial of CD-Chex Plus BC. If the values are still outside the expected range, contact Streck Technical Services at 800-843-0912 or technicalservices@streck.com.
- Clumping of the cell suspension indicates instability or deterioration, in which case the reagent should not be used.

INSTRUCTIONS FOR USE

- Follow instrument manufacturer's instructions for instrument alignment and sample analysis.
- Remove a vial of the control from refrigerator and warm to room temperature (18 °C to 30 °C) for 15 minutes before use.
- Mixing Procedure (**mechanical mixing by vortex or rotator is not recommended**):
For a video demonstration, visit www.streck.com/mixing.

- Holding the vial vertically between the palms of the hands, roll the vial back and forth for 20-30 seconds.

- Hold the vial by the ends between the thumb and finger, and mix by gently inverting the vial at least 8-10 times end-over-end until all cells are thoroughly suspended.

- Aliquot immediately after mixing.
- Subsequent analyses during this test period may be performed by inverting the vial 5 times prior to sampling.

Note: Vials stored for an extended period of time may require extra mixing.

- Return control reagent to refrigeration after sampling to ensure maximum open-vial stability.
- Add recommended monoclonal antibodies according to manufacturer's instructions to each tube and mix gently.
- Incubate according to antibody manufacturer's instructions.
- Add recommended amount of RBC lysing agent and follow manufacturer's instructions.
- Analyze by flow cytometry using your laboratory's established protocol.

PROCEDURES

- Instrument Procedure.** Follow instrument manufacturer's instructions for instrument alignment and sample analysis.
- Monoclonal Antibody Procedure.** Use monoclonal antibodies according to manufacturer's instructions for patient samples.
- RBC Lysing Procedure.** Use RBC lysing agent according to manufacturer's instructions for patient samples.

EXPECTED RESULTS

The mean values provided are derived from replicate analyses on properly compensated flow cytometers.^{1,2} The assay values are obtained using reagents recommended by each instrument manufacturer. The expected ranges listed represent estimates of variation due to different laboratories' reagents, instrument performance, maintenance, and operator technique.

It is recommended that an individual laboratory establish its own control means and ranges for each parameter and level analyzed that reflect the laboratory's specific conditions and protocols. These established control means should fall within the published expected ranges. Data collected from interlaboratory quality control programs can be used as a cumulative approach when calculating ranges.

REFERENCES

- Clinical and Laboratory Standards Institute, H43-A2, Clinical flow cytometric analysis of neoplastic hematology cells. Approved Guideline - Second Edition.
- Clinical and Laboratory Standards Institute, H42-A2, Enumeration of immunologically defined cell populations by flow cytometry. Approved Guideline - Second Edition.

QUALITY CONTROL PROGRAM

Streck offers STATS®, an interlaboratory quality control program, to all customers at no charge. For more information, contact the STATS Department at 800-898-9563 or statsdata@streck.com. Additional information can be found at www.streck.com.

ORDERING INFORMATION

Please call our Customer Service Department at 800-228-6090 for assistance. Additional information can be found online at www.streck.com.

GLOSSARY OF HARMONIZED SYMBOLS

EC REP Authorized Representative in the European Community	LOT Batch Code	Biological Risk	REF Catalog Number	Use By
IVD In Vitro Diagnostic Medical Device	Manufacturer	Consult Instructions For Use	Temperature Limitation	
Glossary of symbols may contain symbols not used in the labeling of this product.				

The brand and product names of the instruments are trademarks of their respective holders.

See www.streck.com/patents for patents that may be applicable to this product.

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