

## **FOR REFERENCE PURPOSES**

**This manual is for Reference Purposes Only. DO NOT use this protocol to run your assays. Periodically, optimizations and revisions are made to the kit and protocol, so it is important to always use the protocol included with the kit.**

**NEXTflex™ BRCA1 & BRCA2 Amplicon Panel**  
(Illumina Compatible)  
Catalog #4220-01 (8 reactions)



## **This product is for research use only.**

**This manual is proprietary to Bioo Scientific Corp., and intended only for customer use in connection with the product(s) described herein and for no other purpose. This document and its contents shall not be used or distributed for any other purpose without the prior written consent of Bioo Scientific. Periodic optimizations and revisions are made to kit components and manuals. Follow the protocol included with the kit.**

**Bioo Scientific makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. There is no warranty of merchantability for this product, or of the fitness of the product for any purpose. Bioo Scientific shall not be liable for any damages, including special or consequential damages, or expense arising directly or indirectly from the use of this product.**

**Bioo Scientific, NEXTflex, AIR, The NGS Experts, qRNA, and NanoQ are trademarks or registered trademarks of Bioo Scientific. All other brands and names contained herein are the property of their respective owners.**

# NEXTflex™ BRCA1 & BRCA2 Amplicon Panel - 4220-01

<b>GENERAL INFORMATION</b> .....	<b>2</b>
Product Overview.....	2
Contents, Storage and Shelf Life.....	2
Required Materials not Provided.....	3
Warnings and Precautions.....	4
<b>NEXTflex™ BRCA1 &amp; BRCA2 AMPLICON PANEL PREPARATION</b> .....	<b>5</b>
NEXTflex™ BRCA1 & BRCA2 Amplicon Panel Preparation Flow Chart.....	5
Starting Material.....	5
Reagent Preparation.....	5
STEP A: PCR I – <i>BRCA1</i> and <i>BRCA2</i> Amplification.....	6
STEP B: PCR I Cleanup.....	7
STEP C: PCR II Amplification (Indexing Reaction) and Cleanup.....	8
<b>LIBRARY VALIDATION</b> .....	<b>10</b>
<b>APPENDIX</b> .....	<b>11</b>
Oligonucleotide Sequences.....	11
Reverse Primer Index Sequences and Reverse Complements.....	11
Low Level Multiplexing.....	11
<b>RELATED PRODUCTS</b> .....	<b>12</b>

## Product Overview

The NEXTflex™ BRCA1 & BRCA2 Amplicon Panel produces Illumina™ compatible, barcoded amplicon libraries in two PCR steps. Libraries are constructed using high quality genomic DNA extracted from blood or cell samples. FFPE or cfDNA samples are not compatible with this kit. This panel contains a total of 130 primer pairs in two pools that allow for the amplification and sequencing of all *BRCA1* and *BRCA2* coding exons. Amplicon regions of interest range between 101 – 229 bp. The regions of interest plus primer pad sites, which comprise the read portion of the libraries, range between 152 – 284 bp. These target regions are amplified and barcoded in two PCR steps (Fig. 1). NEXTflex™ BRCA Amplicon Primer Mixes have been optimized to achieve high coverage uniformity and reduction of off-target reads.

The NEXTflex™ BRCA1 & BRCA2 Amplicon Panel covers 20.4 kilobases comprising 48 coding *BRCA1* and *BRCA2* exons. Libraries have 100% uniformity at  $\geq 0.2x$  mean coverage and > 99% on-target reads. Up to 384 samples can be multiplexed with at least 100x coverage on a single Illumina 2x300 MiSeq lane for detection of germline mutations. Up to 14 samples can be multiplexed at 2000x coverage for detection of somatic mutations. Standard Illumina sequencing primers may be used with this kit.

## Contents, Storage and Shelf Life

The NEXTflex™ BRCA1 & BRCA2 Amplicon Panel contains enough material to prepare 8 samples for Illumina® compatible sequencing. The shelf life of all reagents is 12 months when stored properly. All components should be stored at -20°C, except the Nuclease-free Water and Resuspension Buffer, which can be safely stored at room temperature, and NEXTflex™ Cleanup Beads, which should be stored at 4°C.

Kit Contents	Amount
<b>ORANGE CAP</b>	
NEXTflex™ BRCA Amplicon Primer Mix 1	32 $\mu$ L
<b>RED CAP</b>	
NEXTflex™ BRCA Amplicon Primer Mix 2	32 $\mu$ L
<b>YELLOW CAP</b>	
NEXTflex™ PCR II Barcoded Primer Mix 1 – 8	4 $\mu$ L each
<b>GREEN CAP</b>	
NEXTflex™ PCR Master Mix	288 $\mu$ L
<b>WHITE CAP</b>	
Nuclease-free Water	1.5 mL
Resuspension Buffer	1 mL
<b>BROWN CAP</b>	
NEXTflex™ Cleanup Beads	(2) 900 $\mu$ L

## Required Materials not Provided

- 20-100 ng of high-quality genomic DNA (two 10-50 ng aliquots in up to 34  $\mu$ L nuclease-free water each)
- 96 well PCR Plate Non-skirted (Phenix Research, Cat # MPS-499) or similar
- Adhesive PCR Plate Seal (BioRad, Cat # MSB1001)
- Magnetic Stand -96 (Ambion, Cat # AM10027) or similar
- Thermocycler
- 2, 10, 20, and 200  $\mu$ L pipettes / multichannel pipettes
- Nuclease-free barrier pipette tips
- Vortex
- 80% Ethanol, freshly prepared (room temperature)

## Revision History

Version	Date	Description
v15.11	November 2015	Initial Product Launch
v16.02	February 2016	NEXTflex™ Cleanup Beads are now included with this kit.
v16.04	April 2016	PCR I Reaction and Cleanup volumes have been optimized. PCR I Reactions 1 & 2 are now pooled in PCR II Reaction.

## Warnings and Precautions

Bioo Scientific strongly recommends that you read the following warnings and precautions. Periodically, optimizations and revisions are made to the components and manual. Therefore, it is important to follow the protocol included with the kit. If you need further assistance, you may contact your local distributor or Bioo Scientific at [nextgen@biooscientific.com](mailto:nextgen@biooscientific.com).

- Do not use the kit past the expiration date.
- Ensure pipettes are properly calibrated as PCR is highly sensitive to pipetting error.
- Try to maintain a laboratory temperature of 20°–25°C (68°–77°F).
- Genomic DNA sample quality may vary between preparations. It is the user's responsibility to utilize high quality genomic DNA. Genomic DNA that is heavily nicked or damaged may cause library preparation failure. Absorbance measurements at 260 nm are commonly used to quantify DNA and 260 nm / 280 nm ratios of 1.8-2.0 usually indicate relatively pure DNA. Other quantification methods using fluorescent dyes may also be used. The user should be aware that contaminating RNA, nucleotides and single-stranded DNA may affect the amount of usable DNA in a sample preparation.

## NEXTflex™ BRCA1 & BRCA2 Amplicon Panel Preparation Flow Chart




Figure 1: Sample flow chart with the approximate time necessary for each step.

### Starting Material

The NEXTflex™ BRCA1 & BRCA2 Amplicon Panel has been optimized and validated using 20-100 ng of high-quality genomic DNA (Fig. 2). This kit is not compatible with FFPE or cfDNA samples.

### Reagent Preparation

1. Thaw all tubes on ice.
2. Briefly mix and centrifuge each component just prior to use to ensure material has not lodged in the cap or side of tube.
3. Allow NEXTflex™ Cleanup Beads to come to room temperature and vortex the beads until liquid appears homogenous before every use.

# STEP A: PCR I - *BRCA1* and *BRCA2* Amplification

## Materials

### Bioo Scientific Supplied

**ORANGE CAP** - NEXTflex™ BRCA Amplicon Primer Mix 1

**RED CAP** - NEXTflex™ BRCA Amplicon Primer Mix 2

**GREEN CAP** - NEXTflex™ PCR Master Mix

**WHITE CAP** - Nuclease-free Water

### User Supplied

Thermocycler

96 Well PCR Plate

Adhesive PCR Plate Seal

Ice

For each reaction, 10-50 ng of genomic DNA in up to 34  $\mu$ L water

1. For each sample prepare two separate reactions in adjacent wells of a 96-well PCR Plate on ice as described below. Note: It is recommended to combine these reagents as a master mix if processing multiple samples.

#### Reaction 1

_ $\mu$ L	Genomic DNA (10-50 ng)
_ $\mu$ L	Nuclease-free Water
4 $\mu$ L	NEXTflex™ BRCA Amplicon Primer Mix 1
12 $\mu$ L	NEXTflex™ PCR Master Mix
<hr/>	
50 $\mu$ L	TOTAL

#### Reaction 2

_ $\mu$ L	Genomic DNA (10-50 ng)
_ $\mu$ L	Nuclease-free Water
4 $\mu$ L	NEXTflex™ BRCA Amplicon Primer Mix 2
12 $\mu$ L	NEXTflex™ PCR Master Mix
<hr/>	
50 $\mu$ L	TOTAL

2. Mix thoroughly by pipette.
3. Apply adhesive PCR plate seal and place in thermocycler for the following PCR cycles:

2 min	98°C	Repeat for a total of 9 cycles
<hr/>		
20 sec	98°C	
4 min	62°C	

4. Proceed to Step B: PCR I Cleanup.



## STEP B: PCR I Cleanup

### Materials

#### Bioo Scientific Supplied

WHITE CAP - Resuspension Buffer

BROWN CAP - NEXTflex™ Cleanup Beads (room temperature)

#### User Supplied

80% Ethanol, freshly prepared (room temperature)

96 Well PCR Plate

Magnetic Stand

PCR I Reaction 1 & 2 (from Step A)

1. Add 75  $\mu\text{L}$  of NEXTflex™ Cleanup Beads to each reaction. Mix thoroughly until homogenized.
2. Incubate at room temperature for 5 minutes.
3. Place the 96 well PCR plate on the magnetic stand at room temperature for 5 minutes, or until the supernatant for each sample appears completely clear.
4. Remove and discard the supernatant. Do not disturb beads. Some liquid may remain in wells.
5. With plate on stand, gently add 200  $\mu\text{L}$  of freshly prepared 80% ethanol to each magnetic bead pellet and incubate plate at room temperature for 30 seconds. Carefully remove ethanol by pipette. Repeat this step for a total of 2 ethanol washes. Ensure all ethanol has been removed.
6. Remove the plate from the magnetic stand and let dry at room temperature for 5 minutes or until bead pellets are visibly dry.
7. Resuspend dried beads with 20  $\mu\text{L}$  of Resuspension Buffer. Mix thoroughly by pipette. Ensure beads are no longer attached to the side of the well.
8. Incubate resuspended beads at room temperature for 3 minutes.
9. Place the 96 well PCR plate on the magnetic stand at room temperature for 5 minutes, or until samples appear clear.
10. Transfer 18  $\mu\text{L}$  of clear supernatant from each sample to a new well.
11. Proceed to Step C: PCR II Amplification (Indexing) and Cleanup.

# STEP C: PCR II Amplification (Indexing) and Cleanup

## Materials

### Bioo Scientific Supplied

**YELLOW CAP** - NEXTflex™ PCR II Barcoded Primer Mix (1-8)

**GREEN CAP** - NEXTflex™ PCR Master Mix

**WHITE CAP** - Resuspension Buffer

**BROWN CAP** - NEXTflex™ Cleanup Beads (room temperature)

### User Supplied

Thermocycler

96 Well PCR Plate

Adhesive PCR Plate Seal

80% Ethanol, freshly prepared (room temperature)

Magnetic Stand

Ice

Purified PCR I Reactions 1 & 2 (from step B)

1. For each sample, combine the following reagents on ice in the 96-well PCR plate:

18 µL	Purified PCR I Reaction 1
18 µL	Purified PCR I Reaction 2
2 µL	NEXTflex™ PCR II Barcoded Primer Mix
12 µL	NEXTflex™ PCR Master Mix
<hr/>	
50 µL	TOTAL


2. Mix thoroughly by pipette.
3. Apply adhesive PCR plate seal and place in thermocycler for the following PCR cycles:

2 min	98°C	
<hr/>		
30 sec	98°C	
30 sec	60°C	Repeat for a total of 21 cycles
30 sec	72°C	
<hr/>		
4 min	72°C	


4. Remove plate from thermocycler. Add 40 µL of NEXTflex™ Cleanup Beads to each sample. Mix thoroughly until homogenized.
5. Incubate at room temperature for 5 minutes.
6. Place the 96 well PCR plate on the magnetic stand at room temperature for 5 minutes, or until the supernatant appears completely clear.
7. Remove and discard the supernatant. Do not disturb beads. Some liquid may remain in wells.

8. With plate on stand, gently add 200  $\mu\text{L}$  of freshly prepared 80% ethanol to each magnetic bead pellet and incubate plate at room temperature for 30 seconds. Carefully remove ethanol by pipette. Repeat this step for a total of 2 ethanol washes. Ensure all ethanol has been removed.
9. Remove the plate from the magnetic stand and let dry at room temperature for 5 minutes or until bead pellet is visibly dry.
10. Resuspend dried beads with 22  $\mu\text{L}$  of Resuspension Buffer. Mix thoroughly by pipette. Ensure beads are no longer attached to the side of the well.
11. Incubate resuspended beads at room temperature for 3 minutes.
12. Place the 96 well PCR plate on the magnetic stand at room temperature for 5 minutes, or until sample appears clear.
13. Transfer 20  $\mu\text{L}$  of clear supernatant to a new well.
14. To ensure optimal cluster generation, library concentration must be determined. To determine the library dilution factor required for template preparation, Qubit and/or Bioanalyzer analyses are necessary (Fig 2).

# LIBRARY VALIDATION



**Figure 2. NEXtflex™ BRCA 1 & BRCA 2 Amplicon Panel Libraries.** Libraries were visualized using Agilent Bioanalyzer. High Sensitivity DNA Chip Ladder / Electropherogram:  
A) 100 ng input NEXtflex™ BRCA1 & BRCA2 Amplicon Panel Library (Bioanalyzer gel image)  
B) 100 ng input NEXtflex™ BRCA1 & BRCA2 Amplicon Panel Library (electropherogram)



**Figure 3. Performance of 130 amplicons from NEXtflex™ BRCA 1 & BRCA 2 Amplicon Panel for detection of germline mutations.**

## Oligonucleotide Sequences

NEXTflex™ PCR II Barcoded Primer Mix	
NEXTflex™	Sequence 5' → 3'
PCR II Forward Primer	AATGATACGGCGACCACCGAGATCTACACTCTTTCCCTACACGACGCTCTTCCGATCT
PCR II Reverse Barcoded Primer	CAAGCAGAAGACGGCATAACGAGATXXXXXXXXXXXXXX'GTGACTGGAGTT CAGACGTGTGCTCTTCCGATCT

'XXXXXXXXXXXXXX' denotes the index region of primer. The index sequences and the respective reverse complement sequences contained in each primer are listed below. The reverse complement is the sequence reported in the index read.

## Reverse Primer Index Sequences and Reverse Complements

Barcoded Primer	Sequence 5' → 3'	Reverse Complement
1	GGCCGGCTAGAT	ATCTAGCCGGCC
2	AAGGAAGAGATA	TATCTTTCCTT
3	GGACGGCATCTA	TAGATGCCGTCC
4	AAGGAAGGAGCG	CGCTCCTTCTT
5	GGACGGCGCTCG	CGAGCGCCGTCC
6	CCGACTTCTCGA	TCGAGAGTCCGG
7	GGCCGGCCGAGC	GCTCGGCCGGCC
8	CCGACTGAGCT	AGCTCAGTCCGG

## Low Level Multiplexing

Every combination of sequential odd and even numbered barcodes is fully color balanced at all positions of the index. For example, barcodes 5 and 6 offer opposite colors at every position, but barcodes 6 and 7 do not. Larger pools can be made by combining multiple sets of color balanced pairs. For pools of odd numbers of samples, any barcode can be added to a balanced pool. For example, for a pool of 3 samples, pooling barcodes 5, 6, and any other barcode is acceptable.

For a complete electronic list of the BED and FASTA files for this kit, please follow the instructions on the label on the inside of the kit box.

## RELATED PRODUCTS

### Illumina Compatible RNA NGS Kits and Adapters

NEXTflex™ Rapid Directional RNA-Seq Kit

NEXTflex™ RNA-Seq Barcodes

NEXTflex-96™ RNA-Seq Barcodes

NEXTflex™ Rapid Directional qRNA-Seq™ Kit

NEXTflex™ Small RNA Sequencing Kit v2

NEXTflex™ Small RNA Barcode Primers

NEXTflex™ Poly(A) Beads

### Illumina Compatible DNA NGS Kits and Adapters

NEXTflex™ 16S V4 Amplicon-Seq Kit

NEXTflex™ 16S V4 Amplicon-Seq Kit 2.0

NEXTflex™ 16S V1-V3 Amplicon-Seq Kit

NEXTflex™ 18S ITS Amplicon-Seq Kit

NEXTflex™ Rapid DNA-Seq Kit

NEXTflex™ Cell Free DNA-Seq Kit

NEXTflex™ DNA Barcodes

NEXTflex-96™ DNA Barcodes

NEXTflex-HT™ Barcodes

NEXTflex™ Dual-Indexed DNA Barcodes

NEXTflex™ Bisulfite-Seq Kit

NEXTflex™ Bisulfite-Seq Barcodes

NEXTflex™ Methyl-Seq 1 Kit

NEXTflex™ Msp 1

NEXTflex™ ChIP-Seq Kit

NEXTflex™ ChIP-Seq Barcodes

NEXTflex-96™ ChIP-Seq Barcodes

NEXTflex™ Pre-Capture Combo Kit

NEXTflex™ Rapid Pre-Capture Combo Kit

NEXTflex™ DNA Barcode Blockers

NEXTflex™ PCR-Free DNA Sequencing Kit

NEXTflex™ PCR-Free Barcodes



## **WE WANT TO HEAR FROM YOU!**

Your feedback is important to us. Tell us what you think of our kits by scanning the QR code or visiting our website at [www.biooscientific.com/NGSfeedback](http://www.biooscientific.com/NGSfeedback).

We can't wait to hear from you!



**BIOO SCIENTIFIC®**

**THE NGS EXPERTS™**

Bioo Scientific Corporation · 7050 Burleson Road, Austin, Texas 78744 · [BiooScientific.com](http://BiooScientific.com)  
P: 1.888.208.2246 · F: 512.707.8122 · Bioo Research Products Group · [nextgen@biooscientific.com](mailto:nextgen@biooscientific.com)  
Made in the USA