

Asuragen®
AmplideX®
 PCR/CE *SMN1/2* Kit*




Figure 2. Clinical samples tested with the AmplideX PCR/CE *SMN1/2* Kit have 100% agreement with previous characterization methods.




Figure 3. Detection and differentiation of *SMN1* and *SMN2* hybrid peaks.

KIT ORDERING INFORMATION	
AmplideX PCR/CE <i>SMN1/2</i> Kit*	50 reactions [P/N: A00001]

KIT COMPONENTS	
2X PCR Mix	<i>SMN</i> Control
<i>SMN1/2</i> HEX Primer Mix	Diluent
<i>SMN</i> Calibrator	ROX 1000 Size

93 Blood Samples



Figure 4. Excellent concordance of results between AmplideX and other methods in well-characterized clinical samples.

* For Research Use Only. Not for use in Diagnostic procedures.

† CE-IVD for US Export Only.

^ Two Samples exhibited *SMN2*-to-*SMN1* conversions that were not detected by the AmplideX PCR/CE *SMN1* Kit, consistent with kit limitations.