


Immunophenotyping Control



CD-Chex Select® is the first and only commercially available flow cytometry control with the selected group of CD markers used for immunophenotyping. It is an assayed positive procedural control for reagent staining, erythrocyte lysis, sample preparation, and instrument performance. CD-Chex Select complements the full line of Streck CD-Chex flow cytometry controls.

Convenience

CD-Chex Select eliminates the need for labs to hold and validate patient samples to use as a control for the evaluated markers. Saves lab time and money.

Assayed Instruments

BD Biosciences and Beckman Coulter® flow cytometry systems.

Features and Benefits

- First and only commercially available control with selected group of CD Markers
 - **Sample 1** provides % positive values for CD26, CD41, CD57, CD64, FMC7, TCR α/β, TCR γ/δ.
 - **Sample 2** provides % positive values for CD58, CD61, intracellular CD79a, and intracellular MPO.

Sample 1 and 2 are not sold separately.
- Assayed positive procedural control that can be handled like a patient sample
- Established means and ranges provide the user confidence in instrument results

Stability (Store 2 °C - 10 °C)

Closed-vial stability	90 days
Open-vial stability	30 days

Packaging/Kit Configuration

CD-Chex Select is available in 3.0 mL plastic cap-pierceable vials. Compatible with automated sample preparation instruments.

STATS®

Interlaboratory Quality Control program features online accessibility with enhanced real-time reporting and improved account management capability. Open to all customers at no charge, the Streck STATS program includes personalized, easy-to-read reports that detail performance, identify trends and facilitate real-time prompt corrective action. More information can be found at www.streck.com.

Ordering Information

Description	Catalog #
2 x 3.0 mL (Sample 1, Sample 2)	213554
4 x 3.0 mL (Sample 1, Sample 2)	213555

800.843.0912 | STRECK.COM

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES.

320598-2 2017-07

The brand and product names of the instruments are trademarks of their respective holders.



www.rafer.es

- | | | | | | | | | |
|--|---|--|---|--|---|---|---|---|
| Barcelona
93 645 50 28
barcelona@rafer.es | Bilbao
94 499 85 80
bilbao@rafer.es | La Coruña
981 93 89 26
galicia@rafer.es | Madrid
91 365 15 70
madrid@rafer.es | Málaga
639 359 792
malaga@rafer.es | Sevilla
954 369 334
sevilla@rafer.es | Valencia
96 340 48 00
levante@rafer.es | Zaragoza
976 23 74 00
rafer@rafer.es | Lisboa
21 154 19 98
lisboa@rafer.es |
|--|---|--|---|--|---|---|---|---|