

# MAMMATYPER<sup>®</sup>


Innovación para diagnóstico en cáncer mama

Test de diagnóstico molecular para la detección cuantitativa (qPCR) de la expresión mRNA de ERBB2, ESR1, PGR y MKI67 en cáncer de mama


## MammaTyper Test Result

ERBB2	ESR1	PGR	MKI67	Subtype
Positive	Negative	Negative	Positive	HER2 positive (non-luminal) Subtype according to St Gallen classification 2013*


Details Test Result




ERBB2 Expression		
40-ddCq value	Positive	Cut-off
42.75	3.1 x above Cut-off	41.10



ESR1 Expression		
40-ddCq value	Negative	Cut-off
35.46	5.8 x below Cut-off	38.00



PGR Expression		
40-ddCq value	Negative	Cut-off
33.53	3.9 x below Cut-off	35.50



MKI67 Expression		
40-ddCq value	Positive	Cut-off
36.56	2.1 x above Cut-off	35.50