

F en cemento

セメント + WPowder

Sampling 5mg

AQF-100

Temp 1000/1100deg.C

Absorbent H₂O₂

Volume 10ml


Ion chromatograph

Dionex DX-120

Column AS-12A

Injection 20ul

Suppressor SRS



Nota!

La medida se realiza añadiendo tungsteno como aditivo.

	Sample size (mg)	Results (%)	Average
Cement	5.44	1.33	1.34
	5.03	1.34	
Cement w/Tungsten	6.24	1.92	1.90
	4.97	1.88	

Estimated Value 2%