



AllType™ FASTplex™ NGS Assay

Combining the convenience of multiplex PCR with workflow simplicity

11 loci multiplex



single tube workflow*



*Early pooling of all samples into a single tube within 1st hour

Ordering Information

Product Name	Tests	Cat. No
<i>For Research Use Only. Not for use in diagnostic procedures.</i>		
AllType™ FASTplex™ NGS 11 Loci Flex Kit - 96 ¹	96	ALL-FAST11LFX
AllType™ FASTplex™ NGS 11 Loci Flex Kit - 48 ²	96	ALL-FAST11LX

¹ Includes 96 unique barcodes and all reagents required for amplification and library preparation with Illumina

² Includes 48 unique barcodes and all reagents required for amplification and library preparation with Ion Torrent



TDX-MKT-0104 | OCTOBER 2020 | REV 03

Find out more at onelambda.com

© 2020 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.



www.rafer.es

Barcelona
93 645 50 28
barcelona@rafer.es

Bilbao
94 499 85 80
bilbao@rafer.es

La Coruña
981 93 89 26
galicia@rafer.es

Madrid
91 365 15 70
madrid@rafer.es

Málaga
639 359 792
malaga@rafer.es

Sevilla
954 369 334
sevilla@rafer.es

Valencia
96 340 48 00
levante@rafer.es

Zaragoza
976 23 74 00
rafer@rafer.es

Lisboa
21 154 19 98
lisboa@rafer.es